## 5-ring conductivity measuring cells

Check the measuring cell for transport damage when unpacking.

### Storage

The measuring cell should be stored dry. It should be rinsed well with distilled water immediately before use.

### Measurement

During the measurement, ensure that all five rings are immersed into the measuring solution and that no air bubbles adhere to the platinum rings. To achieve good mixing, the lower tip of the measuring cell should be several millimeters above the vessel bottom, thus leaving enough space for a magnetic stirring bar, for example. However, during the measurement the stirrer should be switched off.

Each conductivity measurement must be carried out at constant temperature. It is measured with the temperature sensor integrated in the 5-ring conductivity measuring cell.

### Cleaning

The measuring cell should be cleaned with distilled water or a 20 % ethanol solution. Before doing this, unscrew and remove the sleeve. Avoid any scratching of the platinum surface. After the cleaning, screw on the sleeve up to the mechanical stop.

### Calibration

Because of alterations of the surface roughness or contamination of the Pt surface, the cell constant changes with time. It should therefore be redetermined at regular intervals.

To do this, use standard solutions of defined conductivity corresponding to the measuring range (see table 1), e.g. Metrohm 6.2301.060 conductivity standard solution (12.88 mS/cm at 25 °C) and 6.2324.000 (100 µS/cm at 25 °C). The calibration should be carried out at 25 °C, or at a temperature at which the conductivity of the standard solution being used is known. As the conductivity of a salt solution is strongly temperature-dependent, the temperature must be held constant to within ±0.1 °C.

The cell constant can be calculated as follows from the measured conductance:

$$c = \frac{\kappa}{G}$$

$\kappa$ : cell constant [1/cm]  
 $\kappa$ : conductivity [S/cm]  
 $G$ : conductance [S]

Calibration and measurement at very low conductivities (< 20 µS/cm) must be carried out using flow-through cells.

Tab. 1: Standards to be used for different measuring ranges

Meas. range [µS/cm]	Conductivity 25 °C / [µS/cm]	Order no.
1 – 1'000	100	6.2324.000
> 1'000	12'880	6.2301.060

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